

# **Why chaperone vectors?**

**A protein folding initiative**

**An open discussion with structural biologists**

# Protein Structure Initiative: Pilot Phase

- ▶ *“Whether the pilot phase achieved its goal depends on how we measure SUCCESS”* (Editorial, PSI-phase 1 and beyond, Nature Structural & Molecular Biology, 11, p. 201, March 2004)
- ★ **Significant increase in throughput of determination of novel protein structures**
- ★ **PSI structures are dominated by structures of single domains primarily from procaryotic proteins**
- ★ **Proteins which “misbehave” are left aside, which leaves out nearly all important proteins for understanding human diseases and for developing cures**
- ▶ *“...it is still unclear how the bottlenecks for eucaryotic and membrane protein structure determination will be overcome”* (ibid)

# Protein Structure Initiative: Update


- *“The availability of suitable recombinant protein is still a major bottleneck in protein structure analysis.”* (Bussow K. et al., Microb Cell Fact., 4: 21, 2005)
- *“Contrary to the popular assumption, the rate of growth of structural data has slowed, and the Protein Data Bank (PDB) has not been growing exponentially since 1995.”* (Levitt M., Proc Natl Acad Sci USA, 104: 3183-3188, 2007)
- ★ **PDB growth rates are steadily decreasing since 1997.**
- ★ **Proteins which “misbehave” are left aside, which leaves out nearly all important proteins for understanding human diseases and for developing cures**
- ★ **The number of novel structures is growing largely through computation.**
- ★ **Such structures do not provide sufficient resolution to allow drug discovery through “docking” of potential drug candidates**




# Inherent PSI problems

- ★ Many eucaryotic proteins require specific sets of molecular chaperones for their folding
- ★ Chaperones are required at compatible concentrations with the overexpressed proteins, however the recombinant protein synthesis exceeds the host cell protein folding capacity


**Protein overexpression**



**Misfolded protein**



*Inclusion bodies*




*Protein degradation products*

- *Proteins need to be overexpressed, since significant amounts of pure proteins are required for structural studies*
- *An overexpressed protein could be misfolded in any, even a mammalian system, unless the chaperones are co-expressed at comparable with the protein expression level*

# Such problems are not unique to structural biology

- ★ Only ~10% of human drug target proteins produced in the most widely used expression system, *E.coli*, are suitable for drug screening

*Protein Expression Consortium*  
([www.lifesensors.com/alliances/pxconsortium/pdf](http://www.lifesensors.com/alliances/pxconsortium/pdf))



## ***Common problems...***


- ★ Insolubility
- ★ low expression level
- ★ cytotoxicity
- ★ proteolysis
- ★ poor antigenicity
- ★ poor biological activity

 ***...The root cause is target protein misfolding***


# Solution to protein misfolding

## Chaperone Vectors

Overexpression of both a recombinant protein and a specific set of human molecular chaperones



Chaperone-assisted protein folding




*Correctly folded recombinant protein*

# Prior art: Improvement in protein folding using molecular chaperones

Experiments in insects cells co-infected with recombinant baculoviruses expressing human target proteins and human chaperones


## *Effect of Hsp40 and Hsp70 chaperones on folding a target protein*



### **8-fold increase**

in yield of soluble cytoplasmic target protein BZLF (Yokohama et al., Biochim Biophys Acta, 1493: 119-124, 2000).

## *Effect of Calreticulin on folding target glycoproteins*



### **9-fold increase**


in lipoprotein lipase enzymatic activity (Zhang et al., J. Biol. Chem., 278: 29344-51, 2003).

**2-3 fold increase** in production of secreted soluble HLA-DR4 tetramers (Fourneau et al., J. Immunol. Methods, 285: 253-64, 2004).



# Chaperones' effect is enhanced if they are delivered in the same vector with target protein (chaperone vector)


## 2 vector system



## Modest positive effect (2-fold)


Target protein and chaperone are delivered into host cells from different vectors

## Chaperone Vector



## Strong positive effect (8-fold)


Target protein and chaperone are delivered into host cells from the same vector




PDI was reported to have modest positive effect on a target protein folding in a 2 vector system (Hsu et al., Protein Expr. Purif., 7, 281-8, 1996). However, the positive effect becomes more pronounced if both PDI and the target protein are co-expressed using the same recombinant baculovirus (chaperone vector). In our experiments, improvement in secretion of cysteine-rich E2 glycoprotein was only about 2 times when insect cells were co-infected with 2 recombinant baculoviruses-one expressing E2 glycoprotein and another expressing human PDI. However, improvement was much more significant (about 8-fold) when both the target protein and the PDI were expressed from the same recombinant baculovirus (Belyaev A.S., unpublished).

## Prior Art

# Multi-vector system



Co-introduce vectors into host cells



Good protein yield and folding in some of the cells


Poor folding

Poor yield


## AB Vector

U.S. Patent No. 7,226,781

# Chaperone vector



Introduce vector into host cells



Good yield and folding of a target protein in every cell of the population

# Why chaperone vectors are more efficient

## Multi-vector system

### ★ *Unbalanced*

Some host cells provide for good target protein yield and folding

Variable production of chaperons and target protein in individual cells as:

- a) Individual host cells receive variable number of vectors
- b) The vectors which enter the cells earlier are replicating faster than the vectors, which enter the cells later

### ★ *Wasteful*

Host cell resources are diverted for replication of several vectors

### ★ *Complex*

A combination of several vectors are used to express molecular chaperones and target protein

## Chaperone vector

### ★ *Balanced*

All host cells provide for good target protein yield and folding

Guaranteed synthesis of target protein and molecular chaperones at the defined ratio in all the cells of the population

### ★ *Economical*

Host cell resources are more efficiently employed for the synthesis of target protein and chaperones. The waste is minimized as there is only one vector

### ★ *Simple*

One vector is used to express molecular chaperones and target protein

# Chaperone vector design






## Vector selection:

### Baculovirus vectors

- ★ Nearly as powerful as *E.coli*
- ★ Higher eucaryotic system (insect cells)
- ★ Well developed technology

## Chaperone selection:

Major human chaperones with demonstrated capacity for target protein folding in insect cells

-  **Hsp40**
  -  **Hsc70**
  -  **Calreticulin-** glycoprotein folding
  -  **Bip-** non-glycoprotein folding in the ER
  -  **PDI-** correct formation of disulfide bonds
- Major cytoplasmic chaperones, work in concert


## Conventional vectors

Baculovirus DNA






Polyhedrin promoter drives powerful target protein synthesis

## Marker selection:

 **GFP-** convenient operation

## Humanized chaperone vectors


Powerful expression of molecular chaperones to enhance target protein folding



Polyhedrin promoter drives powerful target protein synthesis

# ProFold™-C1 – for protein folding in cytoplasm


## BacPAK6 ProFold-C1



Protein expression profiles obtained using conventional BacPAK6 vector and ProFold™-C1 chaperone vector. SDS-PAGE of cell extracts infected with recombinant baculoviruses. Coomassie blue staining.

# ProFold™-ER1 – vector for glycoproteins

BacPAK6    ProFold-ER1



Target protein

Target protein is expressed as efficiently as using parent BaculoGold-BacPAK6 vectors.

PDI

PDI drives correct formation of disulfide bonds and assists protein folding. Human PDI can improve yield of secreted target proteins from mammalian, insect, yeast and even *E.coli* cells.

Calreticulin


Calreticulin is a major chaperone assisting membrane or secreted glycoprotein folding. It is secreted from insect cells reaching up to 200 mg/l of the protein in the cell culture media, thus driving production of a co-expressed secreted glycoprotein (Fourneau et al. 2004).

GFP

GFP marker for convenient operation.

Protein expression profiles obtained using conventional BacPAK6 vector and ProFold™-ER1 chaperone vector. SDS-PAGE of insect cell extracts infected with recombinant baculoviruses. Coomassie blue staining.

# Improvement in target protein folding using a chaperone vector



Extracts of insects cells infected with recombinant baculoviruses. SDS-PAGE, Coomassie blue staining. A cytoplasmic target protein was expressed using ProFold™-C1 (+ chaperones) or using a conventional vector BacPAK6 (-chaperones). NC-cells infected with negative control recombinant baculovirus, which does not express target proteins, chaperones or markers.

# How chaperone vectors compare to conventional vectors?

- The same powerful expression of target proteins as conventional vectors
- The same convenience of making recombinant clones using standard kits
- The same tags for protein purification and detection are available
- More convenient to operate due to the GFP marker

★ *Improved quality of your target protein and even more convenience*

★ *Any reason not to switch to the chaperone vectors?*



# Concluding remarks

★ **Chaperone vectors are commercially available**  
([www.abvector.com](http://www.abvector.com))

★ **Ready to go, just insert ORF encoding target protein into convenient cloning sites**

★ **Collaborations are welcome with the focus on “misbehaving” novel eucaryotic proteins**



**AB Vector**

*www.abvector.com E-mail: [info@abvector.com](mailto:info@abvector.com) Tel: 1-866-683-2867 Fax: 1-760-454-2465*